

Cresko Laboratory Procedures and Protocols

William A. Cresko

Saturday, December 2, 2023

Table of contents

How to use this book	4
I General Laboratory Protocols	5
1 Contact Information	6
II Molecular Protocols	7
2 cDNA basic	8
2.1 Introduction	8
2.2 Materials:	8
2.3 Solutions:	8
2.4 Procedure:	8
2.4.1 First strand synthesis	8
2.4.2 Reaction can be scaled up to accommodate more starting RNA	9
3 2x Turbo	10
3.1 Introduction	10
3.2 Materials:	10
3.3 Solutions:	10
3.4 Procedure:	11
III Vertebrate Husbandry	12
4 Twenty Gallon Aquarium Cleaning	13
4.1 Introduction	13
4.2 Materials:	13
4.3 Solutions:	13
4.4 Procedure:	14
4.5 Air difuser cleaning:	14

IV Daphnia Husbandry	15
5 Placeholder_Daphnia	16
5.1 xxx	16
5.1.1 xxx	16
V Bioinformatic	17
6 A field guide to base R	19
6.1 Introduction	19
6.1.1 Prerequisites	19
6.2 Selecting multiple elements with [.	20
6.2.1 Subsetting vectors	20
6.3 Summary	21
7 Summary	22
References	23
Appendices	24
A Appendix 1 - Sbf1 Barcodes in 96 Well Plate	24
B Appendix 2	32

How to use this book

This is a Quarto book that contains all of the Procedures and Protocols for the Cresko Laboratory in the Institute of Ecology and Evolution at the University of Oregon.

The book is organized into major sections that contain

- General Laboratory Protocols or the lab
- More detailed Laboratory Protocols
- Husbandry protocols for vertebrate animals primarily stickleback and pipefish, but also zebrafish
- Husbandry protocols for *Daphnia*
- Bioinformatic protocols including how to get on to **Talapas**

You can scroll through the book using the index on the left, but also use the search field to find all relevant protocols.

There are also useful appendices at the end, as well as a section for the references cited throughout the book.

This book was written in Markdown using Quarto. To learn more about Quarto books visit <https://quarto.org/docs/books>.

Part I

General Laboratory Protocols

1 Contact Information

Col1	Col2	Col3
Cresko Laboratory	541-346-5189	Phone
Bill Cresko	541-285-5446	Cell
Mark Currey	541-505-0006	Cell
Susie Bassham	xxxx	Cell

Part II

Molecular Protocols

2 cDNA basic

2.1 Introduction

- **Purpose:** This procedure describes how to synthesis cDNA for use with PCR.
- **Procedure Type:** Molecular
- **Species:** N/A

2.2 Materials:

- 2 µl Oligo d(T)23 VN (50 µM, NEB; anchored-dT primer)*
- X µl up to 5 µg total RNA
- 1 µl 10 mM dNTP
- water
- 2 µl 10x RT buffer (Invitrogen)
- 4 µl 25 mM MgCl₂
- 2 µl 0.1 mM DTT – Invitrogen
- 1 µl RNase inhibitor – e.g., RNaseOUT (Invitrogen)
- 1 µl Superscript III reverse transcriptase (200 u/µl – Invitrogen)

2.3 Solutions:

NONE

2.4 Procedure:

2.4.1 First strand synthesis

Combine:

- 2 µl Oligo d(T)23 VN (50 µM, NEB; anchored-dT primer)*
- X µl up to 5 µg total RNA

- 1 μ l 10 mM dNTP mix
- Water (if necessary) to bring total to 10 μ l

Heat to 65°C for 5 min., then ice

Collect contents at bottom of tube by brief centrifugation.

Add:

- 2 μ l 10x RT buffer (Invitrogen)
- 4 μ l 25 mM MgCl₂
- 2 μ l 0.1 mM DTT – Invitrogen
- 1 μ l RNase inhibitor – e.g., RNaseOUT (Invitrogen)
- 1 μ l Superscript III reverse transcriptase (200 u/ μ l – Invitrogen)

Mix by gentle aspiration

- 25°C for 5 min.

2.4.2 Reaction can be scaled up to accommodate more starting RNA

Synthesis: Incubate at 50°C for 50 min.

Inactivation: 85°C for 5 min. Chill on ice, collect contents to bottom by short spin.

Destroy RNA template: 1 μ l RNase H (2 u/ μ l), incubate at 37°C for 20 min.

Proceed to PCR. Depending on expression level, may be able to use a dilution of cDNA as template – try 1:50 dilution in EB, use 2 μ l as template in a 20 μ l reaction. Don't dilute your entire amount of cDNA, as some products may require a higher concentration of template.

3 2x Turbo

3.1 Introduction

- **Purpose:** This procedure describes how to create 2x Turbo PCR mix.
- **Procedure Type:** Molecular
- **Species:** N/A

3.2 Materials:

- 33,000 μ l npH₂O
- 2000 μ l MgSO₄ (100mM)
- 1600 μ l 1M Tris-HCl (pH 8.6)
- 800 μ l 1M KCl
- 800 μ l 1M (NH₄)₂SO₄
- 800 μ l Triton-X 100 (10%)
- 400 μ l DMSO (100 %)
- 120 μ l dATP (100mM)
- 120 μ l dGTP (100mM)
- 120 μ l dTTP (100mM)
- 120 μ l dCTP (100mM)
- 80 μ l 100mg/ml BSA

Total = 40 ml of buffer

3.3 Solutions:

NONE

3.4 Procedure:

- Mix above reagents together
- Place in 1.5 ml ependorph tubes
- Store at -20C

Part III

Vertebrate Husbandry

4 Twenty Gallon Aquarium Cleaning

4.1 Introduction

- **Purpose:** This procedure describes how to clean 20 gallon glass tanks.
- **Procedure Type:** Husbandry
- **Species:**
 - Threespine stickleback, (*Gasterosteus aculeatus*),
 - Gulf pipefish (*Syngnathus scovelli*)

Schedule for Cleaning

Tank cleaning is to be done ONLY Monday - Friday

4.2 Materials:

- Scrub pad or sponge
- Cart (you may or may not want to use)
- Old clothes (this can be messy)
- Personal protection equipment (Splash proof glasses or face shield).

4.3 Solutions:

- **Bleach solution:** Make a 10% bleach solution in a 2 gallon bucket. Add 4.5 L of water. Add 0.5 L of bleach and gently stir.
- **Sodium thiosulfate:** Make a 3% solution of sodium thiosulfate in a separate 2 gallon bucket. Add 5 L of water (to line) and 150g (marked on dispenser) of sodium thiosulfate. Mix

Note: When using bleach and/or sodium thiosulfate. Eye protection is required. Please use splash proof glasses or a face shield when using bleach and sodium thiosulfate.

4.4 Procedure:

- Complete bleaching and cleaning of tank. This needs to be done to each tank every 2 months.
- Remove fish from tank and put them into a clean tank. Tanks that are emptied of fish need to be cleaned and sterilized before another batch of fish can be introduced.
- Drain the tank and remove it from the rack. Clean air diffuser as instructed below.
 - Clean the tank and all parts thoroughly with a scrub pad, taking care not to damage the silicon water seals on the inside (algae should be left if very gentle rubbing will not remove it).
 - Squirt about 10 – 20 mls of bleach into the tank. Wash the bleach water thoroughly around the inside of the tank by hand using a pad or sponge exposing all inside portions of the tank to bleach.
 - Rinse the tank thoroughly with hot tap water. Rinse the tank with sodium thiosulfate, and then rinse it again with hot water. Put a few thiosulfate crystals into the tank and leave it.
 - Reassemble the tank and put it back on the rack. Fill with system water and allow water to recirculate for about 30 minutes before adding fish. Watch fish for 15min to look for any signs of distress.
 - Using a dry erase marker record date/time on the front of the tank when system water is turned back on.
- Initial the check list that you have completed the tank cleaning.

4.5 Air difuser cleaning:

- Remove dirty air diffusers from tanks and rinse with tap water to remove excess algae and debris.
- Place in 10% bleach solution for 15-30 minutes.
- Rinse the corner filters with hot water for 5 and then place into 3% sodium thiosulfate for 5 minutes.
- Rinse with hot water for 5 minutes.
- When cleaned air diffusers are placed back into aquaria, observe fish for 15 min for signs of distress.

Part IV

Daphnia Husbandry

5 Placeholder_Daphnia

5.1 xxx

XXXX

5.1.1 xxx

XXXXX

Part V

Bioinformatic

See Knuth (1984) for additional discussion of literate programming.

6 A field guide to base R

6.1 Introduction

To finish off the programming section, we're going to give you a quick tour of the most important base R functions that we don't otherwise discuss in the book. These tools are particularly useful as you do more programming and will help you read code you'll encounter in the wild.

This is a good place to remind you that the tidyverse is not the only way to solve data science problems. We teach the tidyverse in this book because tidyverse packages share a common design philosophy, increasing the consistency across functions, and making each new function or package a little easier to learn and use. It's not possible to use the tidyverse without using base R, so we've actually already taught you a **lot** of base R functions: from `library()` to load packages, to `sum()` and `mean()` for numeric summaries, to the factor, date, and POSIXct data types, and of course all the basic operators like `+`, `-`, `/`, `*`, `|`, `&`, and `!`. What we haven't focused on so far is base R workflows, so we will highlight a few of those in this chapter.

After you read this book, you'll learn other approaches to the same problems using base R, `data.table`, and other packages. You'll undoubtedly encounter these other approaches when you start reading R code written by others, particularly if you're using StackOverflow. It's 100% okay to write code that uses a mix of approaches, and don't let anyone tell you otherwise!

In this chapter, we'll focus on four big topics: subsetting with `[]`, subsetting with `[[` and `$`, the apply family of functions, and `for` loops. To finish off, we'll briefly discuss two essential plotting functions.

6.1.1 Prerequisites

This package focuses on base R so doesn't have any real prerequisites, but we'll load the tidyverse in order to explain some of the differences.

```
library(tidyverse)
```

6.2 Selecting multiple elements with `[]`

`[]` is used to extract sub-components from vectors and data frames, and is called like `x[i]` or `x[i, j]`. In this section, we'll introduce you to the power of `[]`, first showing you how you can use it with vectors, then how the same principles extend in a straightforward way to two-dimensional (2d) structures like data frames. We'll then help you cement that knowledge by showing how various dplyr verbs are special cases of `[]`.

6.2.1 Subsetting vectors

There are five main types of things that you can subset a vector with, i.e., that can be the `i` in `x[i]`:

1. **A vector of positive integers.** Subsetting with positive integers keeps the elements at those positions:

```
x <- c("one", "two", "three", "four", "five")
x[c(3, 2, 5)]
```

```
[1] "three" "two"   "five"
```

By repeating a position, you can actually make a longer output than input, making the term “subsetting” a bit of a misnomer.

```
x[c(1, 1, 5, 5, 5, 2)]
```

```
[1] "one"  "one"  "five" "five" "five" "two"
```

2. **A vector of negative integers.** Negative values drop the elements at the specified positions:

```
x[c(-1, -3, -5)]
```

```
[1] "two"  "four"
```

3. **A logical vector.** Subsetting with a logical vector keeps all values corresponding to a TRUE value. This is most often useful in conjunction with the comparison functions.

```
x <- c(10, 3, NA, 5, 8, 1, NA)

# All non-missing values of x
x[!is.na(x)]
```

```
[1] 10 3 5 8 1
```

```
# All even (or missing!) values of x
x[x %% 2 == 0]
```

```
[1] 10 NA 8 NA
```

Unlike `filter()`, NA indices will be included in the output as NAs.

4. **A character vector.** If you have a named vector, you can subset it with a character vector:

```
x <- c(abc = 1, def = 2, xyz = 5)
x[c("xyz", "def")]
```

```
xyz def
5 2
```

As with subsetting with positive integers, you can use a character vector to duplicate individual entries.

5. **Nothing.** The final type of subsetting is nothing, `x[]`, which returns the complete `x`. This is not useful for subsetting vectors, but as we'll see shortly, it is useful when subsetting 2d structures like tibbles.

6.3 Summary

In this chapter, we've shown you a selection of base R functions useful for subsetting and iteration. Compared to approaches discussed elsewhere in the book, these functions tend to have more of a “vector” flavor than a “data frame” flavor because base R functions tend to take individual vectors, rather than a data frame and some column specification. This often makes life easier for programming and so becomes more important as you write more functions and begin to write your own packages.

This chapter concludes the programming section of the book. You've made a solid start on your journey to becoming not just a data scientist who uses R, but a data scientist who can *program* in R. We hope these chapters have sparked your interest in programming and that you're looking forward to learning more outside of this book.

7 Summary

In summary, this book has no content whatsoever.

`1 + 1`

[1] 2

References

Knuth, Donald E. 1984. “Literate Programming.” *Comput. J.* 27 (2): 97–111. <https://doi.org/10.1093/comjnl/27.2.97>.

A Appendix 1 - Sbf1 Barcodes in 96 Well Plate

well	Barcode	Name (top)	Final top sequence	well	Name (bottom)	Final bottom sequence
A1	AAACGG	SbfI-AAACGG-top	ACACTCTTTTCACCTACACCG	A1	SbfI-AAACGG-bot	ATCTTATAGGATCGCA
A2	AACGTT	SbfI-AACGTT-top	ACACTCTTTTCACCTACACCG	A2	SbfI-AACGTT-bot	ATCTTATAGGATCGCA
A3	AACTGA	SbfI-AACTGA-top	ACACTCTTTTCACCTACACCG	A3	SbfI-AACTGA-bot	ATCTTATAGGATCGCA
A4	AAGACG	SbfI-AAGACG-top	ACACTCTTTTCACCTACACCG	A4	SbfI-AAGACG-bot	ATCTTATAGGATCGCA
A5	AAGCTA	SbfI-AAGCTA-top	ACACTCTTTTCACCTACACCG	A5	SbfI-AAGCTA-bot	ATCTTATAGGATCGCA
A6	AATATC	SbfI-AATATC-top	ACACTCTTTTCACCTACACCG	A6	SbfI-AATATC-bot	ATCTTATAGGATCGCA
A7	AATGAG	SbfI-AATGAG-top	ACACTCTTTTCACCTACACCG	A7	SbfI-AATGAG-bot	ATCTTATAGGATCGCA
A8	ACAAGA	SbfI-ACAAGA-top	ACACTCTTTTCACCTACACCG	A8	SbfI-ACAAGA-bot	ATCTTATAGGATCGCA
A9	ACAGCG	SbfI-ACAGCG-top	ACACTCTTTTCACCTACACCG	A9	SbfI-ACAGCG-bot	ATCTTATAGGATCGCA
A10	ACATAC	SbfI-ACATAC-top	ACACTCTTTTCACCTACACCG	A10	SbfI-ACATAC-bot	ATCTTATAGGATCGCA

well	Barcode	Name (top)	Final top sequence	well	Name (bottom)	Final bottom sequence
A11	ACCATG	SbfI- ACCATG- top	ACACTCTTTTCACCTACACCG	A11	SbfI- ACCATG- bot	SMIGCTCT/EP66/ATCTAGGTAGATGGCAA
A12	ACCCCC	SbfI- ACCCCC- top	ACACTCTTTTCACCTACACCG	A12	SbfI- ACCCCC- bot	SMIGCTCT/EP66/AGGTACCCACCAATCCCAA
B1	ACTCTT	SbfI- ACTCTT- top	ACACTCTTTTCBCTACACCG	B1	SbfI- ACTCTT- bot	SMIGCTCT/EP66/ATCTACTAGATGGCAA
B2	ACTGGC	SbfI- ACTGGC- top	ACACTCTTTTCBCTACACCG	B2	SbfI- ACTGGC- bot	SMIGCTCT/EP66/AGCTACTAGATGGCAA
B3	AGCCAT	SbfI- AGCCAT- top	ACACTCTTTTCBCTACACCG	B3	SbfI- AGCCAT- bot	SMIGCTCT/EP66/ATCTAGCTAGATGGCAA
B4	AGCGCA	SbfI- AGCGCA- top	ACACTCTTTTCBCTACACCG	B4	SbfI- AGCGCA- bot	SMIGCTCT/EP66/ATCTAGCTAGATGGCAA
B5	AGGGTC	SbfI- AGGGTC- top	ACACTCTTTTCBCTACACCG	B5	SbfI- AGGGTC- bot	SMIGCTCT/EP66/AGATAGCTAGATGGCAA
B6	AGGTGT	SbfI- AGGTGT- top	ACACTCTTTTCBCTACACCG	B6	SbfI- AGGTGT- bot	SMIGCTCT/EP66/ATCTAGCTAGATGGCAA
B7	AGTAGG	SbfI- AGTAGG- top	ACACTCTTTTCBCTACACCG	B7	SbfI- AGTAGG- bot	SMIGCTCT/EP66/ATCTAGTAGATGGCAA
B8	AGTTAA	SbfI- AGTTAA- top	ACACTCTTTTCBCTACACCG	B8	SbfI- AGTTAA- bot	SMIGCTCT/EP66/ATCTAGTAGATGGCAA
B9	ATAGTA	SbfI- ATAGTA- top	ACACTCTTTTCBCTACACCG	B9	SbfI- ATAGTA- bot	SMIGCTCT/EP66/ATCTTTTAGATGGCAA
B10	ATCAAA	SbfI- ATCAAA- top	ACACTCTTTTCBCTACACCG	B10	SbfI- ATCAAA- bot	SMIGCTCT/EP66/ATCTTTTAAATGGCAA
B11	ATGCAC	SbfI- ATGCAC- top	ACACTCTTTTCBCTACACCG	B11	SbfI- ATGCAC- bot	SMIGCTCT/EP66/ATCTCTAGATGGCAA

well	Barcode	Name (top)	Final top sequence	well	Name (bottom)	Final bottom sequence
B12	ATGTTG	SbfI- ATGTTG- top	ACACTCTTTTCC12TACACCG	B12	SbfI- ATGTTG- bot	ACACTCTTTTCC12TACACCG
C1	ATTCCG	SbfI- ATTCCG- top	ACACTCTTTTCC1TACACCG	C1	SbfI- ATTCCG- bot	ACACTCTTTTCC1TACACCG
C2	CAAAAA	SbfI- CAAAAA- top	ACACTCTTTTCC2TACACCG	C2	SbfI- CAAAAA- bot	ACACTCTTTTCC2TACACCG
C3	CAATCG	SbfI- CAATCG- top	ACACTCTTTTCC3TACACCG	C3	SbfI- CAATCG- bot	ACACTCTTTTCC3TACACCG
C4	CACCTC	SbfI- CACCTC- top	ACACTCTTTTCC4TACACCG	C4	SbfI- CACCTC- bot	ACACTCTTTTCC4TACACCG
C5	CAGGCA	SbfI- CAGGCA- top	ACACTCTTTTCC5TACACCG	C5	SbfI- CAGGCA- bot	ACACTCTTTTCC5TACACCG
C6	CATACT	SbfI- CATACT- top	ACACTCTTTTCC6TACACCG	C6	SbfI- CATACT- bot	ACACTCTTTTCC6TACACCG
C7	CCATTT	SbfI- CCATTT- top	ACACTCTTTTCC7TACACCG	C7	SbfI- CCATTT- bot	ACACTCTTTTCC7TACACCG
C8	CCCGGT	SbfI- CCCGGT- top	ACACTCTTTTCC8TACACCG	C8	SbfI- CCCGGT- bot	ACACTCTTTTCC8TACACCG
C9	CCCTAA	SbfI- CCCTAA- top	ACACTCTTTTCC9TACACCG	C9	SbfI- CCCTAA- bot	ACACTCTTTTCC9TACACCG
C10	CCGAGG	SbfI- CCGAGG- top	ACACTCTTTTCC10TACACCG	C10	SbfI- CCGAGG- bot	ACACTCTTTTCC10TACACCG
C11	CCGCAT	SbfI- CCGCAT- top	ACACTCTTTTCC11TACACCG	C11	SbfI- CCGCAT- bot	ACACTCTTTTCC11TACACCG
C12	CCTAAC	SbfI- CCTAAC- top	ACACTCTTTTCC12TACACCG	C12	SbfI- CCTAAC- bot	ACACTCTTTTCC12TACACCG

well	Barcode	Name (top)	Final top sequence	well	Name (bottom)	Final bottom sequence
D1	CGAGGC	SbfI- CGAGGC- top	ACACTCTTTTCCTACACCG	D1	SbfI- CGAGGC- bot	ACACTCTTTTCCTACACCG
D2	CGCAGA	SbfI- CGCAGA- top	ACACTCTTTTCCTACACCG	D2	SbfI- CGCAGA- bot	ACACTCTTTTCCTACACCG
D3	CGCGTG	SbfI- CGCGTG- top	ACACTCTTTTCCTACACCG	D3	SbfI- CGCGTG- bot	ACACTCTTTTCCTACACCG
D4	CGGTCC	SbfI- CGGTCC- top	ACACTCTTTTCCTACACCG	D4	SbfI- CGGTCC- bot	ACACTCTTTTCCTACACCG
D5	CGTCTA	SbfI- CGTCTA- top	ACACTCTTTTCCTACACCG	D5	SbfI- CGTCTA- bot	ACACTCTTTTCCTACACCG
D6	CGTGAT	SbfI- CGTGAT- top	ACACTCTTTTCCTACACCG	D6	SbfI- CGTGAT- bot	ACACTCTTTTCCTACACCG
D7	CTACAG	SbfI- CTACAG- top	ACACTCTTTTCCTACACCG	D7	SbfI- CTACAG- bot	ACACTCTTTTCCTACACCG
D8	CTCGCC	SbfI- CTCGCC- top	ACACTCTTTTCCTACACCG	D8	SbfI- CTCGCC- bot	ACACTCTTTTCCTACACCG
D9	CTGCGA	SbfI- CTGCGA- top	ACACTCTTTTCCTACACCG	D9	SbfI- CTGCGA- bot	ACACTCTTTTCCTACACCG
D10	CTGGTT	SbfI- CTGGTT- top	ACACTCTTTTCCTACACCG	D10	SbfI- CTGGTT- bot	ACACTCTTTTCCTACACCG
D11	CTTATG	SbfI- CTTATG- top	ACACTCTTTTCCTACACCG	D11	SbfI- CTTATG- bot	ACACTCTTTTCCTACACCG
D12	CTTTGC	SbfI- CTTTGC- top	ACACTCTTTTCCTACACCG	D12	SbfI- CTTTGC- bot	ACACTCTTTTCCTACACCG
E1	GAAATG	SbfI- GAAATG- top	ACACTCTTTTCCTACACCG	E1	SbfI- GAAATG- bot	ACACTCTTTTCCTACACCG

well	Barcode	Name (top)	Final top sequence	well	Name (bottom)	Final bottom sequence
E2	GAACCA	Sbfl- GAACCA- top	ACACTCTTTTCE2TACACCG	E2	SMIGCTCT/EP66/ATGCTTACGATGGCAA	ATGCTTACGATGGCAA
E3	GACGAC	Sbfl- GACGAC- top	ACACTCTTTTCE3TACACCG	E3	SMIGCTCT/EP66/AGTTGTCAGATGGCAA	AGTTGTCAGATGGCAA
E4	GACTCT	Sbfl- GACTCT- top	ACACTCTTTTCE4TACACCG	E4	SMIGCTCT/EP66/ATGAGTCTAGTGGCAA	ATGAGTCTAGTGGCAA
E5	GAGAGA	Sbfl- GAGAGA- top	ACACTCTTTTCE5TACACCG	E5	SMIGCTCT/EP66/ATCTGTACAGATCGCAA	ATCTGTACAGATCGCAA
E6	GATCGT	Sbfl- GATCGT- top	ACACTCTTTTCE6TACACCG	E6	SMIGCTCT/EP66/ATCGATACGATCGCAA	ATCGATACGATCGCAA
E7	GCAGAT	Sbfl- GCAGAT- top	ACACTCTTTTCE7TACACCG	E7	SMIGCTCT/EP66/ATCTGCAGATCGCAA	ATCTGCAGATCGCAA
E8	GCATGG	Sbfl- GCATGG- top	ACACTCTTTTCE8TACACCG	E8	SMIGCTCT/EP66/ATCAGCATGGTCGCAA	ATCAGCATGGTCGCAA
E9	GCCGTA	Sbfl- GCCGTA- top	ACACTCTTTTCE9TACACCG	E9	SMIGCTCT/EP66/ATATGGCATGATGGCAA	ATATGGCATGATGGCAA
E10	GCGACC	Sbfl- GCGACC- top	ACACTCTTTTCE10TACACCG	E10	SMIGCTCT/EP66/AGGTGGGAGATGGCAA	AGGTGGGAGATGGCAA
E11	GCGCTG	Sbfl- GCGCTG- top	ACACTCTTTTCE11TACACCG	E11	SMIGCTCT/EP66/ATCTGGGGCTGTCGCAA	ATCTGGGGCTGTCGCAA
E12	GCTCAA	Sbfl- GCTCAA- top	ACACTCTTTTCE12TACACCG	E12	SMIGCTCT/EP66/ATCTAGTACGATGGCAA	ATCTAGTACGATGGCAA
F1	GGACTT	Sbfl- GGACTT- top	ACACTCTTTTCE1TACACCG	E1	SMIGCTCT/EP66/ATCTTCAGCTATCGCAA	ATCTTCAGCTATCGCAA
F2	GGCAAG	Sbfl- GGCAAG- top	ACACTCTTTTCE2TACACCG	E2	SMIGCTCT/EP66/ATCTGGCAAGTTCGCAA	ATCTGGCAAGTTCGCAA
		bot			bot	

well	Barcode	Name (top)	Final top sequence	well	Name (bottom)	Final bottom sequence
F3	GGGCGC	Sbfl- GGGCGC- top	ACACTCTTTTCC E3 TACACCG	SMIGCTCT/EP66/ ATCTGGCGCAGATCGG AA	bot	
F4	GGGGCG	Sbfl- GGGGCG- top	ACACTCTTTTCC E4 TACACCG	SMIGCTCT/EP66/ ATCTGGCGCAGATCGG AA	bot	
F5	GGTACA	Sbfl- GGTACA- top	ACACTCTTTTCC E5 TACACCG	SMIGCTCT/EP66/ ATGTGCTAGATCGG AA	bot	
F6	GGTTTG	Sbfl- GGTTTG- top	ACACTCTTTTCC E6 TACACCG	SMIGCTCT/EP66/ ATCAAGCTATGTCGG AA	bot	
F7	GTAAGT	Sbfl- GTAAGT- top	ACACTCTTTTCC E7 TACACCG	SMIGCTCT/EP66/ ATCTCAAGATCGG AA	bot	
F8	GTATCC	Sbfl- GTATCC- top	ACACTCTTTTCC E8 TACACCG	SMIGCTCT/EP66/ AGCTCAATAGGTGG AA	bot	
F9	GTCATC	Sbfl- GTCATC- top	ACACTCTTTTCC E9 TACACCG	SMIGCTCT/EP66/ AGATGATCAGATGG AA	bot	
F10	GTGCCT	Sbfl- GTGCCT- top	ACACTCTTTTCC E10 TACACCG	SMIGCTCT/EP66/ ATGGATGAGATCGG AA	bot	
F11	GTGTAA	Sbfl- GTGTAA- top	ACACTCTTTTCC E11 TACACCG	SMIGCTCT/EP66/ ATCAGATCAATCGG AA	bot	
F12	GTTGGA	Sbfl- GTTGGA- top	ACACTCTTTTCC E12 TACACCG	SMIGCTCT/EP66/ ATCTGATGAGATCGG AA	bot	
G1	TAAGCT	Sbfl- TAAGCT- top	ACACTCTTTTCC G1 TACACCG	SMIGCTCT/EP66/ ATGTTAAGATCGG AA	bot	
G2	TAATTC	Sbfl- TAATTC- top	ACACTCTTTTCC G2 TACACCG	SMIGCTCT/EP66/ AGATTAATGATGG AA	bot	
G3	TACACA	Sbfl- TACACA- top	ACACTCTTTTCC G3 TACACCG	SMIGCTCT/EP66/ ATGTCAAGATCGG AA	bot	

well	Barcode	Name (top)	Final top sequence	well	Name (bottom)	Final bottom sequence
G4	TACGGG	SbfI- TACGGG- top	ACACTCTTTTCG4TACACCGSMIGCTCT/EP66/ATCTTGAAGATCGGAA		TACGGG- bot	
G5	TAGTAT	SbfI- TAGTAT- top	ACACTCTTTTCG5TACACCGSMIGCTCT/EP66/ATCTTGAAGATCGGAA		TAGTAT- bot	
G6	TATCAC	SbfI- TATCAC- top	ACACTCTTTTCG6TACACCGSMIGCTCT/EP66/ATCTTGAAGATCGGAA		TATCAC- bot	
G7	TCAAAG	SbfI- TCAAAG- top	ACACTCTTTTCG7TACACCGSMIGCTCT/EP66/ATCTTGAAGATCGGAA		TCAAAG- bot	
G8	TCCTGC	SbfI- TCCTGC- top	ACACTCTTTTCG8TACACCGSMIGCTCT/EP66/ATCTTGAAGATCGGAA		TCCTGC- bot	
G9	TCGATT	SbfI- TCGATT- top	ACACTCTTTTCG9TACACCGSMIGCTCT/EP66/ATCTTGAAGATCGGAA		TCGATT- bot	
G10	TCGCCA	SbfI- TCGCCA- top	ACACTCTTTTCG10TACACCGSMIGCTCT/EP66/ATCTTGAAGATCGGAA		TCGCCA- bot	
G11	TCGGAC	SbfI- TCGGAC- top	ACACTCTTTTCG11TACACCGSMIGCTCT/EP66/ATCTTGAAGATCGGAA		TCGGAC- bot	
G12	TCTCGG	SbfI- TCTCGG- top	ACACTCTTTTCG12TACACCGSMIGCTCT/EP66/ATCTTGAAGATCGGAA		TCTCGG- bot	
H1	TCTTCT	SbfI- TCTTCT- top	ACACTCTTTTCG13TACACCGSMIGCTCT/EP66/ATCTTGAAGATCGGAA		TCTTCT- bot	
H2	TGAACC	SbfI- TGAACC- top	ACACTCTTTTCG14TACACCGSMIGCTCT/EP66/ATCTTGAAGATCGGAA		TGAACC- bot	
H3	TGACAA	SbfI- TGACAA- top	ACACTCTTTTCG15TACACCGSMIGCTCT/EP66/ATCTTGAAGATCGGAA		TGACAA- bot	
H4	TGCCCCG	SbfI- TGCCCCG- top	ACACTCTTTTCG16TACACCGSMIGCTCT/EP66/ATCTTGAAGATCGGAA		TGCCCCG- bot	

well	Barcode	Name (top)	Final top sequence	well	Name (bottom)	Final bottom sequence
H5	TGCTTA	SbfI- TGCTTA- top	ACACTCTTTTCH6TACACCGSMIGCTCT/EP66/ATCATGCAATATGGCAA		TGCTTA- bot	
H6	TGGGGA	SbfI- TGGGGA- top	ACACTCTTTTCH6TACACCGSMIGCTCT/EP66/ATCTTCGAGCAATCGCAA		TGGGGA- bot	
H7	TTATGA	SbfI- TTATGA- top	ACACTCTTTTCH7TACACCGSMIGCTCT/EP66/ATCTTAAAGATGGCAA		TTATGA- bot	
H8	TTCCGT	SbfI- TTCCGT- top	ACACTCTTTTCH8TACACCGSMIGCTCT/EP66/ATCTGTACAGATGGCAA		TTCCGT- bot	
H9	TTCTAG	SbfI- TTCTAG- top	ACACTCTTTTCH9TACACCGSMIGCTCT/EP66/ATCTGTAAGATGGCAA		TTCTAG- bot	
H10	TTGAGC	SbfI- TTGAGC- top	ACACTCTTTTCH10TACACCGSMIGCTCT/EP66/AGCTTCAAGATGGCAA		TTGAGC- bot	
H11	TTTAAT	SbfI- TTTAAT- top	ACACTCTTTTCH11TACACCGSMIGCTCT/EP66/ATCTATAAATGGCAA		TTTAAT- bot	
H12	TTTGTC	SbfI- TTTGTC- top	ACACTCTTTTCH12TACACCGSMIGCTCT/EP66/AGCTATAAGATGGCAA		TTTGTC- bot	

B Appendix 2

Hah Hah